

# ThermoLyn supra flexible

Processing the thermoplastic material

Technical information 7.5.3



## 1 Introduction

This technical information provides support when working with ThermoLyn supra flexible, a thermoplastic material (EVA).

This document is directed to trained prosthetists. It is a prerequisite that the qualified personnel are trained in the handling of the various materials, machines and tools.

This technical information does not claim to be exhaustive. Reading this technical information does not substitute reading the instructions for use for all required products.

### 1.1 Applications


ThermoLyn supra flexible is permanently elastic, it can be wiped clean and moulded thermoplastically.

The material has a range of applications. The following table shows examples of use, depending on material thickness:

Material thickness	Example of use
2 mm	Flexible proximal sections of orthoses
3 mm	Tabs Flexible inner sockets upper limb prosthetics (child fitting)
4 mm	Flexible inner sockets upper and lower limb prosthetics
5 mm	Flexible inner sockets lower limb prosthetics (hip disarticulation area)
6 mm	
9 mm	Flexible inner sockets lower limb prosthetics (plates for vacuum forming in tensioning frame)
12 mm	
15 mm	

### 1.2 Flowchart

The entire process is shown in the following flowchart. All work steps described in this document are highlighted in bold.



## 2 Preparation

The following preparations must be made in order to work effectively:

- Collecting the materials and tools
  - Components and auxiliary devices

- Materials
- Machines, equipment and accessories
- Tools
- Preparatory work

## 2.1 Collecting the materials and tools

The materials and tools used in the photos within this technical information are listed in the tables below. The prosthetist assumes full responsibility for the use of any other materials.

Tools	
Designation	Reference number
Scissors	–
Measuring devices	–
Vacuum generation unit	–
Vacuum tube (with vacuum sealing disc)	755x104*
Heat gun	756E11=*
PTFE Teflon foil, trimmed	759Y27=3
Deburring knife	718H5
Heat resistant gloves	641H13
Clean silicone for smoothing (e.g. silicone sanding cone or pieces of liners)	–

Materials	
Designation	Reference number
ThermoLyn supra flexible	616T112*
Lamination disc	29Y57
Isopropyl alcohol	634A58
Silicone-based parting agent (e.g. BetaSil or silicone bonding agent)	–
Polyethylene adhesive tape	627B4
Neoprene® cellular rubber tape	627B5*

## 2.2 Plaster model requirements

ThermoLyn supra flexible is processed directly on the plaster model. A stockinette is not used, because the ThermoLyn supra flexible would adhere to it.

The plaster model must meet the following requirements:

- Made from 87G2 porous plaster. Stockinettes are not required, because the plaster is permeable to air.
- Completely dry. A vacuum can only be generated with a completely dry plaster model.  
**TIP:** The plaster model can be dried in a fan oven at max. 60° C.
- Surface completely smoothed and free of bubbles. The material is very soft during processing and even small irregularities are transferred.  
**TIP:** Bubbles can be made good with Cellona plaster bandages. Finally polish with a piece of stockinette.

### INFORMATION

87G2 porous plaster is only required when a vacuum is to be created.

## 2.3 Processing information

The relevant information for the processing of ThermoLyn supra flexible is summarised in this section.

### Heating

### INFORMATION

A protective film is located on the material, which has to be removed prior to heating.

- The material can be processed in the temperature range between **80 °C** and **120 °C**.
- **TIP:** The larger the piece of material to be processed, the lower the temperature that should be selected. The material becomes very soft, making it more difficult to process larger pieces.
- The duration of heating depends on the material thickness and the equipment in which it is heated. Lift the material at one corner to test its condition.
- The material can be heated with the following equipment: hotplate, infrared oven, convection oven

- Place the material on a Teflon film during heating.

### **Welding**

- The material fuses with itself as long as it is at processing temperature.
- The material can be welded subsequently.

### **Smoothing**

- The material can be shaped and smoothed using silicone (e.g. silicone sanding cones, pieces of silicone liners or silicone mats).
- **TIP:** Briefly heat the surface with a hot air gun to restore the shine.

### **Cutting/sanding**

- The material can be cut with scissors or a scalpel.
- The material can be sanded at the socket router (e.g. with the following tools: Habermann cleaner, silicone sanding cone, foam roller).

## **3 Procedure**

Three examples are used to demonstrate the processing of ThermoLyn supra flexible:

- Inner socket for a forearm prosthesis
- Flexible proximal section for an orthosis
- Inner socket with tab for hip disarticulation prosthesis

All relevant processing techniques are demonstrated in these examples. Information on heating the material is provided in the "Processing information" section (see Page 3).

### **3.1 Inner socket forearm prosthesis**

#### **INFORMATION**

The following steps should be carried out by two persons.

#### **INFORMATION**

The material is placed on the plaster model from above and welded on the underneath.

- Align the plaster model so that the shortest side faces downwards. The seam is created on this side.

Select the suitable material thickness.

**Ottobock recommends:** 4 mm

Cut the material to size:

- According to the dimensions of the plaster model and an allowance for the seam.

Remove the protective foil.

Clean the material with isopropyl alcohol.

Place the material on a Teflon film and heat it.



Lift the heated material at its 4 corners and drape it over the plaster model.



**INFORMATION: Do not pull on the material, to prevent the creation of different wall thicknesses.**

Shape the material to the plaster model and join the seam on the underside of the plaster model.

Switch on the vacuum unit.

**INFORMATION: Keep the vacuum unit switched on throughout the entire processing.**



Using scissors, cut off any excess material at the seam. This prevents the weight of the excess material from stretching the area around the seam.



Seal the area around the vacuum tube with sealing tape (e.g. made from Vulkollan®).

Allow the material to cool to room temperature.



Using a scalpel, cut off any excess material flush with the seam.

Heat the seam area with a hot air gun until the material can be shaped.





Smooth the seam with a silicone sanding cone.

**TIP:** Reheat the area around the seam after smoothing, to restore the shine.

### Placing the lamination disks

The lamination disks are embedded in the thermoplastic material. They have an M4 interior thread for screwing the thermoplastic material to a stable frame.



Screw a screw into each lamination disk until the thread is completely filled. This prevents the material from entering the thread.

Place the lamination disk in the required position.



Heat the lamination disk and the surrounding material with the hot air gun.




When the material is soft enough, press the lamination disks into the material. The material is pressed through the slotted holes in the lamination disk.



Press the lamination disk into the material until the material that has been pressed through is flush with the surrounding material.

Remove the screw from the lamination disk.




**TIP:** Infilling is only required if the lamination disk is not covered sufficiently.

Cut off a small strip of material.

Heat the strip and the material surrounding the lamination disk.

Fill in the indentation around the lamination disk with the strip.



Reheat the area to fuse the material.



Smooth the area with a clean piece of silicone (e.g. silicone liner).

**TIP:** Reheat the area after smoothing to restore the shine.

Repeat this process for all lamination disks.

### Finishing the prosthetic socket

Trim and sand the socket edge.

Sand away material over dummies (e.g. for electrodes).

Remove the prosthetic socket from the mould and finish it. Finishing the prosthetic socket is not part of this document.

### 3.2 Flexible proximal section orthosis

#### INFORMATION

The following steps should be carried out by two persons.



Mark the flexible proximal section on the plaster model.  
Select the suitable material thickness.

**Ottobock recommends: 2 mm**

Cut the material to size:

- According to the marked contours and an allowance.

Remove the protective foil.

Clean the material with isopropyl alcohol.

Place the material on a Teflon film and heat it.



Lift the heated material at its 4 corners and drape it over the plaster model.



**INFORMATION: Do not pull on the material, to prevent the creation of different wall thicknesses.**

Mould the heated material to the plaster model. A piece of silicone (e.g. silicone liner) can be used for this purpose.

Allow the material to cool to room temperature.






Transfer the marked flexible proximal section to the material.

Mark the positions of the laminating disks.

### Placing the lamination disks


The lamination disks are embedded in the thermoplastic material. They have an M4 interior thread for screwing the thermoplastic material to a stable frame.



Screw a screw into each lamination disk until the thread is completely filled. This prevents the material from entering the thread.

Place the lamination disk in the required position.

Heat the lamination disk and the surrounding material with the hot air gun.



When the material is soft enough, press the lamination disks into the material until the material extrudes slightly from the slotted holes in the lamination disk.

Remove the screw from the lamination disk.



**TIP:** Infilling is only required if the lamination disk is not covered sufficiently.

Cut off a small strip of material.

Heat the strip and the material surrounding the lamination disk.

Fill in the indentation around the lamination disk with the strip.

Reheat the area to fuse the material.



Smooth the area with a clean piece of silicone (e.g. silicone liner).

**TIP:** Reheat the area after smoothing to restore the shine.

Repeat this process for all lamination disks.

### Finishing of the flexible proximal section



Cut the material according to the marked contours.

**If required:** Bevel the edges to create a softer transition to the laminate.

Attach the material to the plaster model with double-sided adhesive tape.

Carry out the lamination process. The lamination process does not form part of this document.

## 3.3 Inner socket hip disarticulation

### Fabrication of the tab

Select the suitable material thickness.

**Ottobock recommends:** 2 mm

Cut the material to size:

- According to the dimensions of the tab and an allowance.

Remove the protective foil.

Clean the material with isopropyl alcohol.

Place the material on a Teflon film and heat it.



Lift the heated material and place it on the plaster model.



**INFORMATION: Do not pull on the material, to prevent the creation of different wall thicknesses.**

Mould the heated material to the plaster model. A piece of silicone (e.g. silicone liner) can be used for this purpose.

Allow the material to cool to room temperature.

Mark the contour of the tab on the material.



Mark the edge to be bevelled.

Bevel the edge at the socket router.

Place the material on the plaster model.



Mask the part of the tab that is to be welded to the prosthetic socket with PE adhesive tape.



Apply a silicone-based parting agent (e.g. BetaSil, Impresil or silicone bonding agent) to the material.





Distribute the parting agent thinly and evenly e.g. with a piece of foam.



Remove the adhesive tape immediately after the parting agent has been applied.  
Leave the parting agent to harden.



Cut a strip of Neoprene® cellular rubber tape to size and place it where the inner socket is to be opened.  
Attach the Neoprene® cellular rubber tape with a stapler.

### Fabricating the inner socket

#### INFORMATION

The following steps should be carried out by two persons.

#### INFORMATION

The material is placed on the plaster model from above and welded on the underneath.

- Align the plaster model so that the shortest side faces downwards. The seam is created on this side.





Select the suitable material thickness.

**Ottobock recommends: 5 mm**

Cut the material to size:

- According to the dimensions of the plaster model and an allowance for the seam.

Remove the protective foil.

Clean the material with isopropyl alcohol.

Place the material on a Teflon film and heat it.

Lift the heated material at its 4 corners and drape it over the plaster model.



**INFORMATION: Do not pull on the material, to prevent the creation of different wall thicknesses.**

Mould the material to the contours of the plaster model and join the seam on the underside of the plaster model.

- Join the seam on the underside of the plaster model.
- Join the material on the side of the residual limb.



Switch on the vacuum unit.

**INFORMATION: Keep the vacuum unit switched on throughout the entire processing.**

Shape the material to avoid seams or folds in the area of the future prosthetic socket as much as possible.



Using scissors, cut off any excess material at the seams. This prevents the weight of the excess material from stretching the areas.

Seal the area around the vacuum tube with sealing tape (e.g. made from Vulkollan®).



Allow the material to cool to room temperature.

Using a scalpel, cut off any excess material flush with the seam.

Heat the seam area with a hot air gun until it becomes malleable.

Smooth the seam with a clean piece of silicone or a silicone sanding cone.

**TIP:** Reheat the area around the seam after smoothing, to restore the shine.

### Placing the lamination disks

The lamination disks are embedded in the thermoplastic material. They have an M4 interior thread for screwing the thermoplastic material to a stable frame.



Screw a screw into each lamination disk until the thread is completely filled. This prevents the material from entering the thread.

Place the lamination disk in the required position.


Heat the lamination disk and the surrounding material with the hot air gun.



When the material is soft enough, press the lamination disks into the material. The material is pressed through the slotted holes in the lamination disk.

Press the lamination disk into the material until the material that has been pressed through is flush with the surrounding material.

Remove the screw from the lamination disk.



**TIP:** Infilling is only required if the lamination disk is not covered sufficiently.

Cut off a small strip of material.

Heat the strip and the material surrounding the lamination disk.

Fill in the indentation around the lamination disk with the strip.

Reheat the area to fuse the material.



Smooth the area with a piece of silicone (e.g. silicone liner).

**TIP:** Reheat the area after smoothing to restore the shine.



Repeat this process for all lamination disks.


### Finishing the prosthetic socket



Cut open the prosthetic socket in the middle above the Neoprene® cellular rubber tape. Do not cut into the tab while doing so.


Trim the socket edge.

Remove the prosthetic socket from the mould and finish it. Finishing the prosthetic socket is not part of this document.





[illegible]



# Kundenservice/Customer Service

## Europe

Otto Bock HealthCare Deutschland GmbH  
Max-Näder-Str. 15 · 37115 Duderstadt · Germany  
T +49 5527 848-3411 · F +49 5527 848-1414  
healthcare@ottobock.de · www.ottobock.de

Otto Bock Healthcare Products GmbH  
Kaiserstraße 39 · 1070 Wien · Austria  
T +43 1 5269548 · F +43 1 5267985  
vertrieb.austria@ottobock.com · www.ottobock.at

Otto Bock Adria Sarajevo D.O.O.  
Omladinskih radnih brigada 5  
71000 Sarajevo · Bosnia-Herzegovina  
T +387 33 766200 · F +387 33 766201  
obadria@bih.net.ba · www.ottobockadria.com.ba

Otto Bock Bulgaria Ltd.  
41 Tzar Boris III' Blvd. · 1612 Sofia · Bulgaria  
T +359 2 80 57 980 · F +359 2 80 57 982  
info@ottobock.bg · www.ottobock.bg

Otto Bock Suisse AG  
Pilatusstrasse 2 · CH-6036 Dierikon  
T +41 41 455 61 71 · F +41 41 455 61 70  
suisse@ottobock.com · www.ottobock.ch

Otto Bock ČR s.r.o.  
Protetická 460 · 33008 Zruč-Senec · Czech Republic  
T +420 377825044 · F +420 377825036  
email@ottobock.cz · www.ottobock.cz

Otto Bock Iberica S.A.  
C/Majada, 1 · 28760 Tres Cantos (Madrid) · Spain  
T +34 91 8063000 · F +34 91 8060415  
info@ottobock.es · www.ottobock.es

Otto Bock France SNC  
4 rue de la Réunion · CS 90011  
91978 Courtaboeuf Cedex · France  
T +33 1 69188830 · F +33 1 69071802  
information@ottobock.fr · www.ottobock.fr

Otto Bock Healthcare plc  
32, Parsonage Road · Englefield Green  
Egham, Surrey TW20 0LD · United Kingdom  
T +44 1784 744900 · F +44 1784 744901  
bockuk@ottobock.com · www.ottobock.co.uk

Otto Bock Hungária Kft.  
Tatai út 74. · 1135 Budapest · Hungary  
T +36 1 4511020 · F +36 1 4511021  
info@ottobock.hu · www.ottobock.hu

Otto Bock Adria d.o.o.  
Dr. Franje Tuđmana 14 · 10431 Sveta Nedelja · Croatia  
T +385 1 3361 544 · F +385 1 3365 986  
ottobockadria@ottobock.hr · www.ottobock.hr

Otto Bock Italia Srl Us  
Via Filippo Turati 5/7 · 40054 Budrio (BO) · Italy  
T +39 051 692-4711 · F +39 051 692-4710  
info.italia@ottobock.com · www.ottobock.it

Otto Bock Benelux B.V.  
Ekkersrijt 1412 · 5692 AK  
Son en Breugel · The Netherlands  
T +31 499 474585 · F +31 499 476250  
info.benelux@ottobock.com · www.ottobock.nl

Industria Ortopédica Otto Bock Unip. Ltda.  
Av. Miguel Bombarda, 21 - 2º Esq.  
1050-161 Lisboa · Portugal  
T +351 21 3535587 · F +351 21 3535590  
ottobockportugal@mail.telepac.pt

Otto Bock Polska Sp. z o. o.  
Ulica Koralowa 3 · 61-029 Poznań · Poland  
T +48 61 6538250 · F +48 61 6538031  
ottobock@ottobock.pl · www.ottobock.pl

Otto Bock Romania srl  
Șos de Centura Chitila · Mogoșoia Nr. 3  
077405 Chitila, Jud. Ilfov · Romania  
T +40 21 4363110 · F +40 21 4363023  
info@ottobock.ro · www.ottobock.ro

OOO Otto Bock Service  
p/o Pultikovo, Business Park „Greenwood”,  
Building 7, 69 km MKAD  
143441 Moscow Region/Krasnogorskiy Rayon  
Russian Federation  
T +7 495 564 8360 · F +7 495 564 8363  
info@ottobock.ru · www.ottobock.ru

Otto Bock Scandinavia AB  
Koppargatan 3 · Box 623 · 60114 Norrköping · Sweden  
T +46 11 280600 · F +46 11 312005  
info@ottobock.se · www.ottobock.se

Otto Bock Slovakia s.r.o.  
Röntgenova 26 · 851 01 Bratislava 5 · Slovak Republic  
T +421 2 32 78 20 70 · F +421 2 32 78 20 89  
info@ottobock.sk · www.ottobock.sk

Otto Bock Sava d.o.o.  
Industrijska bb · 34000 Kragujevac · Republika Srbija  
T +381 34 351 671 · F +381 34 351 671  
info@ottobock.rs · www.ottobock.rs

Otto Bock Ortopedi ve  
Rehabilitasyon Tekniği Ltd. Şti.  
Ali Dursun Bey Caddesi · Latı Lokum Sokak  
Meriç Sitesi B Blok No: 6/1  
34387 Mecidiyeköy-Istanbul · Turkey  
T +90 212 3565040 · F +90 212 3566688  
info@ottobock.com.tr · www.ottobock.com.tr

## Africa

Otto Bock Algérie E.U.R.L.  
32, rue AHCÈNE Outaleb - Coopérative les Mimosas  
Mackle-Ben Aknoun · Alger · DZ Algérie  
T +213 21 913863 · F +213 21 913863  
information@ottobock.fr · www.ottobock.fr

Otto Bock Egypt S.A.E.  
28 Soliman Abaza St. Mohandessein · Giza · Egypt  
T +202 330 24 390 · F +202 330 24 380  
info@ottobock.com.eg · www.ottobock.com.eg

Otto Bock South Africa (Pty) Ltd  
Building 3 Thornhill Office Park · 94 Bekker Road  
Midrand · Johannesburg · South Africa  
T +27 11 312 1255  
info-southafrica@ottobock.co.za  
www.ottobock.co.za

## Americas

Otto Bock Argentina S.A.  
Av. Belgrano 1477 · CP 1093  
Ciudad Autónoma de Buenos Aires · Argentina  
T +54 11 5032-8201 / 5032-8202  
atencionclientes@ottobock.com.ar  
www.ottobock.com.ar

Otto Bock do Brasil Tecnica Ortopédica Ltda.  
Alameda Maria Tereza, 4036, Bairro Dois Córregos  
CEP: 13.278-181, Valinhos-São Paulo · Brasil  
T +55 19 3729 3500 · F +55 19 3269 6061  
ottobock@ottobock.com.br · www.ottobock.com.br

Otto Bock HealthCare Canada  
5470 Harvester Road  
Burlington, Ontario, L7L 5N5, Canada  
T +1 289 288-4848 · F +1 289 288-4837  
infocanada@ottobock.com · www.ottobock.ca

Otto Bock HealthCare Andina Ltda.  
Calle 138 No 53-38 · Bogotá · Colombia  
T +57 1 8619988 · F +57 1 8619977  
info@ottobock.com.co · www.ottobock.com.co

Otto Bock de Mexico S.A. de C.V.  
Prolongación Calle 18 No. 178-A  
Col. San Pedro de los Pinos  
C.P. 01180 México, D.F. · Mexico  
T +52 55 5575 0290 · F +52 55 5575 0234  
info@ottobock.com.mx · www.ottobock.com.mx

Otto Bock HealthCare  
Two Carlson Parkway North, Suite 100  
Minneapolis, MN 55447 · USA  
T +1 763 553 9464 · F +1 763 519 6153  
usa.customerservice@ottobockus.com  
www.ottobockus.com

## Asia/Pacific

Otto Bock Australia Pty. Ltd.  
Suite 1.01, Century Corporate Centre  
62 Norwest Boulevard  
Baulkham Hills NSW 2153 · Australia  
T +61 2 8818 2800 · F +61 2 8814 4500  
healthcare@ottobock.com.au · www.ottobock.com.au

Beijing Otto Bock Orthopaedic Industries Co., Ltd.  
B12E, Universal Business Park  
10 Jiuxianqiao Road, Chao Yang District  
Beijing, 100015, P.R. China  
T +8610 8598 6880 · F +8610 8598 0040  
news-service@ottobock.com.cn  
www.ottobock.com.cn


Otto Bock Asia Pacific Ltd.  
Unit 1004, 10/F, Greenfield Tower, Concordia Plaza  
1 Science Museum Road, Tsim Sha Tsui  
Kowloon, Hong Kong · China  
T +852 2598 9772 · F +852 2598 7886  
info@ottobock.com.hk · www.ottobock.com

Otto Bock HealthCare India Pvt. Ltd.  
20th Floor, Express Towers  
Nariman Point, Mumbai 400 021 · India  
T +91 22 2274 5500 / 5501 / 5502  
information@indiaottobock.com · www.ottobock.in

Otto Bock Japan K. K.  
Yokogawa Building 8F, 4-4-44 Shibaura  
Minato-ku, Tokyo, 108-0023 · Japan  
T +81 3 3798-2111 · F +81 3 3798-2112  
ottobock@ottobock.co.jp · www.ottobock.co.jp

Otto Bock Korea HealthCare Inc.  
4F Agaworld Building · 1357-74, Seocho-dong  
Seocho-ku, 137-070 Seoul · Korea  
T +82 2 577-3831 · F +82 2 577-3828  
info@ottobockkorea.com · www.ottobockkorea.com

Otto Bock South East Asia Co., Ltd.  
1741 Phaholyothin Road  
Kwaeng Chatuchark · Khet Chatuchark  
Bangkok 10900 · Thailand  
T +66 2 930 3030 · F +66 2 930 3311  
obsea@ottobock.co.th · www.ottobock.co.th



Printed by:

Otto Bock HealthCare GmbH  
Max-Näder-Straße 15 · 37115 Duderstadt · Germany  
T +49 5527 848-0 · F +49 5527 72330  
healthcare@ottobock.de · www.ottobock.com

Ottobock has a certified Quality Management System in accordance with ISO 13485.