

ottobock.

ComforTex

Innovative padding materials



Quality for life

Information for specialist dealers



ComforTex innovative padding materials

Article name	ComforTex soft	ComforTex smooth	ComforTex grippy	ComforTex air
Reference number	623F109	623F110	623F112	623F62
Chemical composition, orthosis side	100% polyamide	100% polyamide	100% polyamide	100% polyamide
Compatible with hook-and-loop, orthosis side	Yes	Yes	Yes	Yes
Chemical composition, skin side	16% elastane 84% polyamide	80% polyamide 20% elastane	84% polyamide 16% elastane	100% polyester
Compatible with hook-and-loop, skin side	No	No	No	No
Area of application	For padding applications in orthotics	For padding applications in orthotics	For padding applications in orthotics	For padding applications in orthotics
Features	<ul style="list-style-type: none"> • Especially soft surface • Excellent padding characteristics • High resilience • Excellent wearer comfort • Individual shaping possible • HF and ultrasonic welding possible • 30°C/86°F gentle cycle 	<ul style="list-style-type: none"> • Especially smooth surface • Excellent padding characteristics • High resilience • Excellent wearer comfort • Individual shaping possible • HF and ultrasonic welding possible • 30°C/86°F gentle cycle 	<ul style="list-style-type: none"> • Special knitting technology prevents slipping in the nap-direction • Excellent padding characteristics • High resilience • Excellent wearer comfort • Individual shaping possible • HF and ultrasonic welding possible • 30°C/86°F gentle cycle 	<ul style="list-style-type: none"> • Breathable • Good padding characteristics • High resilience • High wearer comfort • Individual shaping possible • 30°C/86°F gentle cycle

• ComforTex soft
623F109

• ComforTex smooth
623F110

• ComforTex grippy
623F112

• ComforTex air
623F62

Order information

623F109 ComforTex soft, 623F110 ComforTex smooth, 623F112 ComforTex grippy, 623F62 ComforTex air

Article number	623F*=1-7	623F*=2-7	623F*=5-7
Length	1 m	2 m	5 m
Width	1400 mm	1400 mm	1400 mm
Thickness	4 mm	4 mm	4 mm
Colour	Black/black (7)	Black/black (7)	Black/black (7)

* applies to all ComforTex materials

