

QuickFit Strap QuickFit Buckle

Better orthosis fit and comfort in next to no time


The combination of the QuickFit Strap and QuickFit Buckle is ideal for individually adjustable orthoses and their fabrication. It offers technological solutions that result in an orthosis with a better fit and improved comfort, as it is easy to individually adjust to a user's respective needs.

1 QuickFit Strap

Minute adjustments using only one hand are possible thanks to the BOA closure. Users can therefore influence the compression or hold of the orthosis individually. The lightweight, handy strap also offers a high level of comfort.

2 QuickFit Buckle

The magnetic component ensures effortless locking of the closure. Easy and intuitive opening is possible with one hand.



Benefits at a glance

- Closure and adjustment mechanism for orthoses
- The length can be individually shortened
- The magnetic component ensures that the mechanism is aligned and locked effortlessly
- Minute adjustments can be made on the go
- Easy, one-handed operation by turning the reel
- Quick opening and closing
- Robust and durable
- Lightweight material
- Various combinations of the QuickFit Strap and QuickFit Buckle are possible

Technical data

514Z20 QuickFit Strap

Article number	Length	Width	Colour
514Z20=25-7		approx. 25 mm	
514Z20=38-7	approx. 73–63 cm	approx. 38 mm	black
514Z20=50-7		approx. 50 mm	

514Z21 QuickFit Buckle

Article number	Length	Width	Colour
514Z21=25-7		approx. 25 mm	
514Z21=38-7	approx. 18 cm	approx. 38 mm	black
514Z21=50-7		approx. 50 mm	

BOA® is a registered trademark of Boa Technology, Inc.