



ottobock.

Accessories for socket routers

Quality for life

An Ottobock Group company

jos america

Grinding tools



749F6 Sanding drum

749F6

EquipmentWith conical attachment, 1 Sanding sleeve (Grit 80)



749Y8=* Sanding sleeve

749Y8=73X200X24 749Y8=73X200X40 749Y8=73X200X60

For 749F6 Ottobock Sanding drum

Grit 24 40 60

749Y8=73X200X80 749Y8=73X200X100 749Y8=73X200X150

For 749F6 Ottobock Sanding drum

Grit 80 100 150





749Y16=* Sanding sleeve

	749Y16=A45	749Y16=B45	749Y16=C45
--	------------	------------	------------

For 749Z7=* Ottobock Polisher

Outside diameter Ø 45 mm

Grit 40 80 120

	749Y16=A60	749Y16=B60	749Y16=C60
--	------------	------------	------------

For 749Z7=* Ottobock Polisher

Outside diameter Ø 60 mm

Grit 40 80 120

749Z3=* Habermann Polisher fine, small

	749Z3=5/8	749Z3=M16
Thread	5/8"	M16
Outside diameter Ø	45 mm	



749Y7=* Sanding belt

	749Y7=A45X58	749Y7=B45X58	749Y7=C45X58
For	749Z3=* Habermann Polisher fine, small		
Width	50 mm		
Grit	60	100	150
Outside diameter Ø	45 mm		





749Z2=* Habermann Polisher fine, normal

	749Z2=5/8	749Z2=M16
Thread	5/8"	M16
Outside diameter Ø	60 mm	



749Y7=* Sanding belt

	749Y7=A60X58	749Y7=B60X58	749Y7=C60X58
For	749Z2=* Habermann Polisher fine, normal		
Grit	60	100	150
Outside diameter Ø	60 mm		

749Z4=* Habermann Polisher fine, long

	749Z4=5/8	749Z4=M16
Equipment	With a sanding belt (Grit 100)	
Thread	5/8"	M16
Outside diameter Ø	60 mm	



749Y7=* Sanding belt

	749Y7=A60X100	749Y7=B60X100	749Y7=C60X100
For	749Z4=* Habermann fine polisher, long		
Grit	60	100	150
Outside diameter Ø	60 mm		





749Z5=* Habermann Polisher fine, oversize

	749Z5=5/8	749Z5=M16
Equipment	With a sanding belt (Grit 100)	
Thread	5/8"	M16
Outside diameter Ø	85 mm	



749Y7=* Sanding belt

	749Y7=A85X58	749Y7=B85X58	749Y7=C85X58
For	749Z5=* Habermann Polisher fine, oversize		
Grit	60	100	150
Outside diameter Ø	85 mm		

749Z8=* Sanding drum

	749Z8=5/8"X25	749Z8=M16X25	749Z8=1/2"X25
Thread	5/8"	M16	1/2"-13



749Y25=25 Rubber roll

	749Y25=25
For	749Z8=* Sanding drum



749Y22=* Sanding sleeve

	749Y22=A25	749Y22=C25
For	749Z8=* Sanding drum	
Grit	40 (Coarse)	120 (Fine)





749Z9=* Sanding cone

	749Z9=5/8"	749Z9=M16	749Z9=1/2"
Equipment	With rubber roll, with one sanding sleeve (Grit 80)		
Thread	5/8"	M16	1/2"-13



749Y26=* Sanding sleeve

	749Y26=40	749Y26=80	749Y26=150
For	749Z9 Sanding cone		
Grit	40 (Coarse)	80 (Medium)	150 (Fine)

749Z6=* Sanding drum

	749Z6=5/8"X22	749Z6=M16X30	749Z6=M16X22
Thread	5/8"	M16	M16
Outside diameter Ø	22 mm	30 mm	22 mm
Equipment	With rubber roll, with one sanding sleeve (Grit 40)		

	749Z6=5/8"X30	749Z6=1/2"X22	749Z6=1/2"X30
Thread	5/8"	1/2"-13	1/2"-13
Outside diameter Ø	30 mm	22 mm	30 mm
Equipment	With rubber roll, with one sanding sleeve (Grit 40)		



749Y9=* Sanding sleeve

	749Y9=A22	749Y9=A30	749Y9=C22	749Y9=C30
For	749Z6=* Sanding drum			
Grit	40 (Coarse)		120 (Fine)	
Outside diameter Ø	22 mm	30 mm	22 mm	30 mm

749Y90=40 Sanding sleeve for carbon and compound materials

	749Y90=40
For	749Z6=* Sanding drum
Width	30 mm
Ø	30 mm
Grit	40 (Coarse)

749Z16=* Sanding drum

	749Z16=5/8	749Z16=M16
Thread	5/8"	M16
Equipment	With rubber roll, with one sanding sleeve (Grit 80)	



749Y45=* Sanding sleeve

	749Y45=80	749Y45=150
For	749Z16=* Sanding drum	
Grit	80	150





749Z17=* Sanding drum

	749Z17=5/8	749Z17=M16
Thread	5/8"	M16
Equipment	With rubber roll, with one sanding sleeve (Grit 80)	



749Y48=* Sanding sleeve

	749Y48=80	749Y48=150
For	749Z17=* Sanding drum	
Grit	80	150

749Z18=* Sanding drum

	749Z18=5/8	749Z18=M16
Thread	5/8"	M16
Equipment	With rubber roll, with one sanding sleeve (Grit 80)	



749Y49=* Sanding sleeve

	749Y49=80	749Y49=150
For	749Z18=* Sanding drum	
Grit	80	150





749F8=* Sanding cone

	749F8=5/8X45	749F8=M16X65	749F8=M16X45
Thread	5/8"	M16	M16
Outside diameter Ø	45 mm	65 mm	45 mm
Equipment	With one sanding sleeve (Grit 40)		



	749F8=5/8X65	749F8=1/2"X45	749F8=1/2"X65
Thread	5/8"	1/2"-13	1/2"-13
Outside diameter Ø	65 mm	45 mm	65 mm
Equipment	With one sanding sleeve (Grit 40)		

749Y10=* Sanding sleeve

	749Y10=A45	749Y10=A65	749Y10=B45	749Y10=B65
For	749F8 Sanding cone			
Grit	40		80	
Outside diameter Ø	45 mm	65 mm	45 mm	65 mm
Scope of delivery	6 Pieces			





749F16=* Silicone sanding cones

	749F16=5/8	749F16=M16
Thread	5/8"	M16
Scope of delivery	A set consists of 4 pieces	
To be used for	Polishing and rounding edges of thermoplastic materials	

Showroom Ottobock | Jos America
Industriepark 'het Hoog' · Mandenmaker 14
5253 RC Nieuwkuijk / Netherlands · T +31 (0)73 511 9123
info@josamerica.com · www.ottobock-josamerica.com/en/

Otto Bock HealthCare GmbH
Max-Näder-Str. 15 · 37115 Duderstadt/Germany
peservice@ottobock.com · www.ottobock.com