

Hygiene station

to protect employees and patients

With the new Ottobock hygiene station, you help protect yourself and others from infectious diseases.

The station contains everything necessary for disinfecting the hands and for cleaning and disinfecting orthopaedic devices.

When developing it, special emphasis was placed on interrupting the chain of transmission. The disinfectant and soap dispensers and the tap can be operated with the forearm. In addition, all surfaces are stainless steel so they can be easily cleaned and disinfected.



All benefits at a glance:

- Everything in one place
- Can be used for disinfecting hands and devices
- Smooth surfaces are easy to clean

758E11=1 Hygienestation

Article number	758E11=1
Dimensions W x D x H	1.200 x 643 x 1.785 mm
Equipment	<ul style="list-style-type: none"> • Stainless steel surfaces for easy cleaning • sink with forearm tap • disinfectant and soap dispenser • paper holder • holder for glove boxes • place for the hygiene plan • shelf • instant water heater as an option
Material	Stainless steel (work surface and backsplash), painted steel (base cabinet)
Intended use	Hand disinfection and disinfection of orthopaedic devices



Further information is available at:
www.pem.ottobock.com