

Ottobock Collection Thermopapers and Socket Decors

A piece of individuality



Quality for life



Be inspired!

Thermopapers and socket decors make it possible to enhance surfaces with new designs. They are suitable for a wide range of applications in orthopaedic technology (prosthetics, orthotics). The socket decors permit custom colouring with the use of the lamination technique. Our thermopapers are suitable for giving plastics such as ThermoLyn PE 200, ThermoLyn PP-C and ThermoLyn PP-H an attractive, designer look.

Ottobock offers a varied collection that includes different patterns and colours. No additional workshop equipment or special technical knowledge is required to use them.

In this product information, you will find patterns, colours and design variations that can help you create a device that your patient not only likes to wear but will also want to show off.

Processing instructions



Warm the synthetic material until it becomes completely transparent.



Cut the thermopaper to size and apply it, with the design facing down, onto the pre-cut component; carefully smooth the thermopaper from the middle towards the outside edges.



After approx. 20 seconds, pull the thermopaper off the material.

Benefits at a glance

- ▶ Ottobock thermopapers: for custom colouring of plastics such as ThermoLyn PE 200, ThermoLyn PP-C and ThermoLyn PP-H
- ▶ Environmentally friendly, tear-proof paper
- ▶ Ottobock socket decors: for custom colouring with the lamination technique
- ▶ Proven tubular fabric
- ▶ Large selection of designs to visually alter prostheses and orthoses to meet the patient's preferences
- ▶ No additional workshop equipment or special technical knowledge required

Please note that the colours shown in this product information may differ in actual effect.

Available information materials

Article number	Description
646D696=EN	Product Information for upgrading
646M70=ALL_INT	Design sample folder for thermopapers and socket decors
623P200=1	Thermopaper sample roll
623S200=1	Socket decor sample roll



Pull the decorative fabric around the socket model with the design facing inwards. Using your hand to hold the decorative fabric taut, trace along the contour of the socket using a pencil. Do not use permanent markers, ball-point pens, etc., because these could cause discolouration during the laminating process.




Take in the decorative fabric, leaving 5 mm between the stitches and the inside of the marked line. Allow the distal seam to taper off. This means the seam can be seen only on the medial side. Using scissors, trim the decorative fabric along the seam to a maximum of 2 mm.



Turn the decorative fabric inside out so that the design is facing outwards, and then pull it over the socket model. Place the seam in an inconspicuous position on the socket (e.g. medially).

Ottobock collection of thermopapers



Thermopaper “Tabaluga”


Order example: **article number = length – colour** (623P110=2=5)

Article number	623P110=2-*	623P110=5-*
Length	2 m	5 m
Width	75 cm	75 cm
Colour	Red (2), green (3), blue (5), white (6)	



Thermopaper “Giraffes, light green”

Article number	623P50=2	623P50=5	623P50=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Light green	Light green	Light green



Thermopaper “Starry Sky, blue”

Article number	623P52=2	623P52=5	623P52=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Light blue	Light blue	Light blue



Thermopaper “Bears, red blue yellow”

Article number	623P97=2	623P97=5	623P97=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Red, blue and yellow	Red, blue and yellow	Red, blue and yellow



Thermopaper “On the Road”

Article number	623P98=2	623P98=5	623P98=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Green	Green	Green


Thermopaper “Butterflies”

Article number	623P101=2	623P101=5	623P101=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Purple	Purple	Purple



Thermopaper “Stars, pink”

Article number	623P60=2	623P60=5	623P60=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Pink	Pink	Pink



Thermopaper “Bears, green purple orange”

Article number	623P99=2	623P99=5	623P99=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Green, purple and orange	Green, purple and orange	Green, purple and orange



Thermopaper “Fox”


Article number	623P106=2	623P106=5	623P106=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm



Thermopaper “Bears, blue green grey”


Article number	623P100=2	623P100=5	623P100=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Blue, green and grey	Blue, green and grey	Blue, green and grey






Thermopaper “Kaleidoscope, black red”

Article number	623P92=2	623P92=5	623P92=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Black and red	Black and red	Black and red



Thermopaper “Kaleidoscope, black white”

Article number	623P93=2	623P93=5	623P93=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Black and white	Black and white	Black and white



Thermopaper “Mandala, to colour it yourself”

Article number	623P103=2	623P103=5	623P103=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	For colouring	For colouring	For colouring



Thermopaper “Muscles, white black”

Article number	623P94=2	623P94=5	623P94=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	White and black	White and black	White and black



Thermopaper “Faces”

Article number	623P85=2	623P85=5	623P85=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm


Thermopaper “Raindrops”

Article number	623P87=2	623P87=5	623P87=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Blue	Blue	Blue



Thermopaper “Concrete”

Article number	623P88=2	623P88=5	623P88=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Grey	Grey	Grey



Thermopaper “Camouflage Earth”

Article number	623P86=2	623P86=5	623P86=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Green	Green	Green



Thermopaper “Mini Jeans”

Article number	623P62=2	623P62=5	623P62=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Blue	Blue	Blue



Thermopaper “Mega Jeans”


Article number	623P61=2	623P61=5	623P61=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Blue	Blue	Blue





Thermopaper “Modern Jeans”

Article number	623P84=2	623P84=5	623P84=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Blue	Blue	Blue



Thermopaper “Carbon, sand”

Article number	623P107=2	623P107=5	623P107=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Sand	Sand	Sand



Thermopaper “Carbon, raspberry red”

Article number	623P89=2	623P89=5	623P89=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Raspberry red	Raspberry red	Raspberry red



Thermopaper “Carbon, water blue”

Article number	623P90=2	623P90=5	623P90=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Water blue	Water blue	Water blue



Thermopaper “Carbon, heavy grey”

Article number	623P74=2	623P74=5	623P74=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Dark grey	Dark grey	Dark grey


Thermopaper “Uni Skin Colour”

Article number	623P96=2	623P96=5	623P96=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Beige	Beige	Beige




Thermopaper “Uni Black”

Article number	623P95=2	623P95=5	623P95=10
Length	2 m	5 m	10 m
Width	75 cm	75 cm	75 cm
Colour	Black	Black	Black



Ottobock collection of socket decors



Socket decor “Tabaluga”


Order example: **article number = width – colour** (623S110=15-3)

Article number	623S110=15-*	623S110=20-*
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Red (2), green (3), blue (5), white (6)	



Socket decor “Giraffes, light green”

Article number	623S50=15	623S50=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Light green	



Socket decor “Starry sky, blue”

Article number	623S52=15	623S52=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Light blue	



Socket decor “Bears, red/blue/yellow”

Article number	623S97=15	623S97=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Red, blue and yellow	



Socket decor “On the road”

Article number	623S98=15	623S98=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Green	

Socket decor “Butterflies”

Article number	623S101=15	623S101=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Purple	Purple



Socket decor “Stars, pink”

Article number	623S60=15	623S60=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Pink	Pink



Socket decor “Bears: green, purple, orange”

Article number	623S99=15	623S99=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Green, purple and orange	Green, purple and orange



Socket decor “Fox”


Article number	623S106=15	623S106=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm



Socket decor “Bears, blue, green, grey”


Article number	623S100=15	623S100=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Blue, green and grey	Blue, green and grey






Socket decor “Kaleidoscope, black, red”

Article number	623S92=15	623S92=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Black and red	Black and red



Socket decor “Kaleidoscope, black, white”

Article number	623S93=15	623S93=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Black and white	Black and white



Socket decor “Mandala, colour it yourself”

Article number	623S103=15	623S103=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	For colouring	For colouring



Socket decor “Muscles, white, black”

Article number	623S94=15	623S94=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	White and black	White and black



Socket decor “Faces”

Article number	623S85=15	623S85=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm


Socket decor “Raindrops”

Article number	623S87=15	623S87=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Blue	Blue




Socket decor “Concrete”

Article number	623S88=15	623S88=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Grey	Grey



Socket decor “Camouflage Earth”

Article number	623S86=15	623S86=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Green	Green



Socket decor “Mini jeans”

Article number	623S62=15	623S62=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Blue	Blue



Socket decor “Mega jeans”

Article number	623S61=15	623S61=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Blue	Blue





Socket decor “Modern jeans”

Article number	623S84=15	623S84=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Blue	Blue




Socket decor “Carbon, sand”

Article number	623S107=15	623S107=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Sand	Sand



Socket decor “Carbon, raspberry red”

Article number	623S89=15	623S89=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Raspberry red	Raspberry red



Socket decor “Carbon, water blue”

Article number	623S90=15	623S90=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Water blue	Water blue



Socket decor “Carbon, heavy grey”

Article number	623S74=15-7	623S74=20-7
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Dark grey	Dark grey


Socket decor “Uni skin colour”

Article number	623S96=15	623S96=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Beige	Beige




Socket decor “Uni black”

Article number	623S95=15	623S95=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Black	Black




Socket decor “Ottobock, light blue”

Article number	623S81=15	623S81=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Light blue	Light blue



Socket decor “Ottobock, green”

Article number	623S82=15	623S82=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Green	Green



Socket decor “Ottobock, yellow”

Article number	623S83=15	623S83=20
Length	1.5 m	1.5 m
Width	15 cm	20 cm
Colour	Yellow	Yellow

