

# 749F16 Silicone Sanding Cone Set



Ottobock silicone sanding cones set a new standard for finishing prosthetic sockets and orthoses made of thermoplastic materials. They make polishing and rounding of the edges very simple and clean, resulting in ultimate comfort for the user.

A silicone sanding cone set includes 4 different shapes for the best possible use in the work process:

- 1 Conical with narrow head
- 2 Conical
- 3 Round and cylindrical
- 4 Round and cylindrical with grooves

<b>Reference number</b>	<b>749F16=</b>	<b>5/8</b>	<b>M16</b>
Thread		5/8"	M16